


## NEW PRODUCT

CLIP APPLYING FORCEPS FOR OPEN SURGERY  
SUITABLE FOR ETHICON, PILLING, SLS AND HORIZON CLIPS



Clip applying forceps							
clip size	length of the instrument	ring colour	for Ethicon-Clips	ring colour	for Pilling-Clips	for SLS-Clips	for Horizon-Clips
micro	150 mm			white		<b>204-200</b>	
small	150 mm	blue	<b>204-001</b>	yellow	<b>204-101</b>	<b>204-201</b>	<b>204-301</b>
small-medium	150 mm			purple		<b>204-202</b>	
medium	150 mm	white	<b>204-003</b>	blue	<b>204-103</b>	<b>204-203</b>	<b>204-303</b>
micro	190 mm			white		<b>204-210</b>	
small	190 mm	blue	<b>204-011</b>	yellow	<b>204-111</b>	<b>204-211</b>	<b>204-311</b>
small-medium	190 mm			purple		<b>204-212</b>	
medium	190 mm	white	<b>204-013</b>	blue	<b>204-113</b>	<b>204-213</b>	<b>204-313</b>
medium-large	190 mm	green	<b>204-014</b>	green	<b>204-114</b>	<b>204-214</b>	<b>204-314</b>
large	190 mm	orange	<b>204-015</b>	orange	<b>204-115</b>	<b>204-215</b>	<b>204-315</b>
micro	280 mm			white		<b>204-220</b>	
small	280 mm			yellow		<b>204-221</b>	
small-medium	280 mm			purple		<b>204-222</b>	
medium	280 mm	white	<b>204-023</b>	blue	<b>204-123</b>	<b>204-223</b>	<b>204-323</b>
medium-large	280 mm	green	<b>204-024</b>	green	<b>204-124</b>	<b>204-224</b>	<b>204-324</b>
large	280 mm	orange	<b>204-025</b>	orange	<b>204-125</b>	<b>204-225</b>	<b>204-325</b>